




Kee Diagnostics
because *life* matters

ANA-Chem Touch



- ➔ **Compact & Easy to Operate**
- ➔ **7 inch color-LCD Touch screen**
- ➔ **High Speed Micro CPU**
- ➔ **Real time curve showing**
- ➔ **Profile Based Programming**
- ➔ **Wavelength : 7 Standard Filters, 340, 405, 492, 510, 546
578, 630nm, 1 Free Position**
- ➔ **Memory for 10,000 sample results**
- ➔ **Light source hibernate function making lamp long life time**
- ➔ **10 Incubators (Optional), $37 \pm 0.3^{\circ}\text{C}$**
- ➔ **QC: Levy Jennings Plots**
- ➔ **RS232 interface, PC connecting**



Kee Diagnostics
because life matters

ANA-Chem Touch

PROFILES

Liver Function:	ALT, AST, ALP, T-BIL, D-BIL, TP, ALB, Y-GT
Renal Function:	BUN, CRE, UA, URE, CO2-CP
Lipids:	T-CHO, TG, HDL-C, LDL-C
Cardiac function:	CK, CK-MB, LDH
Immun proteins:	IGA, IGG, IGM, C3
Others:	GLU, HGB
Test mode:	Kinetic, End point, Twopoint, Absorbance
Wavelength:	7 standard filters, 340, 405, 492, 510, 546, 578, 630nm, 1 free position
Wavelength precision:	±2nm, width< 10nm
Absorbance range:	0~4.500 Abs
Absorbance precision:	outside 0.0001 Abs, inside 0.00001 Abs
Carryover:	< 1%
Repeatability CV:	< 1%
Incubator:	37°C, 30°C, 25°C, and others, Precision: ±0.1°C
Sample volume:	0~3000ul(500ul recommend)
Data Storage:	10000"
Quality control curve:	X SD CV% L-J quality control curve for each test item
Flow cell:	32 ul quartz glass, 10mm
Lamp:	6V 10W Halogen Lamp with longer life time
Screen:	7 inch color LCD
CPU:	ARM series high speed micro CPU
Printer:	Internal thermal sensitive printer
Interface:	RS-232
Work Environment:	Temperature 0°C~40°C, Humidity: 80%
Power supply:	100~240VAC, 50~60Hz, 100VA
Dimension:	L340mm* W270mm*H160mm
Weight:	5KG



Kee Diagnostics
because life matters

Kee Diagnostics Pvt. Ltd.
A-8, 3rd Floor, Naraina Industrial Area, Phase II, New Delhi-110028
E Mail Id: info@keediagnosics.in Website: www.keediagnosics.in
Ph: 91 11 43136000