

URIT-670

Microplate Washer

| Specifications |

Plate type: U, V or flat-bottom

Wash head: 8-probe or 12-probe

Wash programs: Up to 50

Wash cycles: 1 ~ 20

Volume range: 50 ~ 450 μ L/well

Residual volume: \leq 2 μ L/well

Dispense precision: CV \leq 3%

Soak Time: 1 ~ 90sec tunable

Absorbing Time: 1-10sec tunable

Shaking Time: 0-60sec tunable

Display: Large LCD display

Keypad: 8 keys

